

## Xtend Cover Protection for Prosthetics

*Xtend Cover is a cosmetic and protective cover for knee-prosthesis. It is designed to give your legs a personal touch so that you can show them with pride. The design gives you a personal expression and becomes an accessory, while also protecting you and your prosthetic leg.*



The protective cover consists of two parts, a front cover and a back cover that are easily connected with two screws in the bottom and a snap-cover on the top. It is based on a 3D printed base, with light and durable polyamide, with a tailored leather-cover.



Rubber spacers keep the rest of the cover in place and is shock-absorbing if you fall or rest on your knees.

### Technical data

Tube-size	30/34 mm
Weight	350 gr
Core material	Polyamide
External material	Natural leather

### Fitting knee-joints

Fits C-leg, Genium, RheoKnee and most of the mechanical knee-joints on the market.

## Ordering information Design Honeycomb



White Leather/Matte Black    Black Leather/Glossy Black    Brown Leather/Silver Metallic    Black Leather/Silver Metallic

Front	Tube size	Order number:
White leather	30 mm	C100-WH-01-30-01
White leather	34 mm	C100-WH-01-34-01
Black leather	30 mm	C100-BL-01-30-01
Black leather	34 mm	C100-BL-01-34-01
Brown leather	30 mm	C100-BR-01-30-01
Brown leather	34 mm	C100-BR-01-34-01



Back	Leg:	Tube size	Ordering number:
Matte black	Right	30 mm	C100-MB-01-30-2L
Matte black	Left	30 mm	C100-MB-01-30-2L
Matte black	Right	34 mm	C100-MB-01-34-2L
Matte black	Left	34 mm	C100-MB-01-34-2L
Glossy black	Right	30 mm	C100-GB-01-30-2L
Glossy black	Left	30 mm	C100-GB-01-30-2L
Glossy black	Right	34 mm	C100-GB-01-34-2L
Glossy black	Left	34 mm	C100-GB-01-34-2L
Silver Metallic	Right	30 mm	C100-SM-01-30-2L
Silver Metallic	Left	30 mm	C100-SM-01-30-2L
Silver Metallic	Right	34 mm	C100-SM-01-34-2L
Silver Metallic	Left	34 mm	C100-SM-01-34-2L



# Ordering information

## Design Flower



White Leather/Matte Black

Black Leather/Glossy Black

Brown Leather/Silver Metallic

Front	Tube size	Ordering number
White leather	30 mm	C100-WH-02-30-01
White leather	34 mm	C100-WH-02-34-01
Black leather	30 mm	C100-BL-02-30-01
Black leather	34 mm	C100-BL-02-34-01
Brown leather	30 mm	C100-BR-02-30-01
Brown leather	34 mm	C100-BR-02-34-01



Back	Leg	Tube size	Order number
Matte black	Right	30 mm	C100-MB-02-30-2R
Matte black	Left	30 mm	C100-MB-02-30-2L
Matte black	Right	34 mm	C100-MB-02-34-2R
Matte black	Left	34 mm	C100-MB-02-34-2L
Glossy black	Right	30 mm	C100-GB-02-30-2R
Glossy black	Left	30 mm	C100-GB-02-30-2L
Glossy black	Right	34 mm	C100-GB-02-34-2R
Glossy black	Left	34 mm	C100-GB-02-34-2L
Silver Metallic	Right	30 mm	C100-SM-02-30-2R
Silver Metallic	Left	30 mm	C100-SM-02-30-2L
Silver Metallic	Right	34 mm	C100-SM-02-34-2R
Silver Metallic	Left	34 mm	C100-SM-02-34-2L



