

Xtend® Foot Först H100

Technical Data Specification

Xtend Foot Först is a patented foot prosthesis, which with its 17 degrees of flexibility, functions almost like a natural foot. The defining feature of Xtend Foot Först is its superior lateral mobility, providing secure compensation on uneven surfaces. A low center of gravity assures unrivaled stability. This is thanks to its unique mix of the adaptable and flexible glass fiber and carbon fiber, a material known for its ability to create energy return.

Recommended L-codes: L-5981 and L-5986



Weight

Lightweight construction.
Ex: Size P426 = 399 grams / 0.88 lbs.



Mounting

Universal pyramid adapter that fits all standard components on the market

Water proof

Guaranteed in fresh-, pool- and salt water.



Material

Glass fiber, carbon fiber and titanium.

Flexible

17° inversion
17° eversion
– comparable with a human foot.



Warranty

Xtend Foot Först – 36 months.
Foot shell – 9 months.



Max weight patient

150 kg / 330 lbs (size category P7).



One unit


Xtend Foot Först has no loose parts.



Pre-assembled by delivery

Xtend Foot Först, Footshell, Spectra sock already in place.

Technical facts Xtend® Foot Först H100



A Foot shell height in mm

B Build height in mm

Xtend Foot Selection tool



Find the perfect fit with our digital Foot Selection tool!

		Foot size (cm)													
		23		24		25		26		27		28		29	
Category		A	B	A	B	A	B	A	B	A	B	A	B	A	B
P3	Max 60 kg/130 lbs	72	106	70	107										
P4	Max 80 kg/175 lbs	72	106	70	107	76	115	82	117	84	118	86	120	87	120
P5	Max 100 kg/220 lbs	72	108	70	108	76	116	82	118	84	118	86	122	87	122
P6	Max 125 kg/275 lbs					76	117	82	118	84	118	86	123	87	123
P7	Max 150 kg/330 lbs							82	117	84	118	86	123	87	124

Size guide

Category		Foot size (cm)						
		23	24	25	26	27	28	29
P3	Max 60 kg/130 lbs	H100-323	H100-324					
P4	Max 80 kg/175 lbs	H100-423	H100-424	H100-425	H100-426	H100-427	H100-428	H100-429
P5	Max 100 kg/220 lbs	H100-523	H100-524	H100-525	H100-526	H100-527	H100-528	H100-529
P6	Max 125 kg/275 lbs			H100-625	H100-626	H100-627	H100-628	H100-629
P7	Max 150 kg/330 lbs				H100-726	H100-727	H100-728	H100-729

How to order

Select the footshell color Caucasian (FC), Brown (FB), Tan (FT) or Jet Black (FJ), and add the code as a suffix. Example: H100-426-FB.

For left (L) or right (R), add L or R as the last part of the suffix. Example H100-426-FBL